

Product Code: 750546

SPILLCHAMP 50 LITRE CHEMICAL SPILL KIT

Convenient Kit for Small to Medium Chemical Spills



OVERVIEW

The SpillChamp 50 Litre Chemical Spill Kit is your reliable solution for managing hazardous chemical spills swiftly and securely. Housed in a durable plastic box, this kit is designed for convenient storage and quick access, allowing for immediate response to minimise risks and cleanup time.

Perfectly suited for handling aggressive spills, it provides peace of mind when working with harsh chemicals like acids. With safety as a priority, the kit includes essential personal protective equipment (PPE) to protect users.

Each kit absorbs approximately 50 litres and includes -

- 25 x chemical absorbent pads
- 6 x chemical absorbent sock
- 1 x PVC waste bag and tie
- 1 x PPE kit (including goggles and gloves)
- 1 x black/yellow barrier tape

PRODUCT FEATURES

- Robust plastic box convenient reusable box to store the kit
- Small to medium capacity suitable for localised chemical spills
- Waste bags and PPE included safe and easy disposal of waste items
- Refillable reusable plastic box, refill with spill contents for repeat use
- Suitable for aggressive chemical, acid or solvent spills in labs, warehouses and manufacturing facilities

TECHNICAL DATA





SPECIFICATION	CHEMICAL SPILL PADS	CHEMICAL SPILL SOCKS
Material:	Hydrophobic textile	
Color:	Yellow	
Item absorbency:	1litre	3 - 4 litres
Quantity:	25	6
Length:	500mm	1200mm
Width:	400mm	70mm
Total kit absorbency:	50 litres	
Liquids absorbed:	Acids and aggressive chemicals	
Container type:	Plastic box	
Applications	Small to medium capacity suitable for localised aggressive chemical, acid or solvent spills in labs, warehouses and manufacturing facilities	



INSTRUCTIONS FOR USE

1. Acquire the spill kit: Bring the kit to the spill location - it is recommended to have the kit close by to places proned to spills

2. Cordon spill area: Prevent access to the area by using barriers or warnings to deter pedestrian entry

3. Apply PPE: Ensure any cleanup staff wear gloves and goggles to avoid injury

4. Apply materials to spill: Apply the pads and socks to the spill and allow to soak as long as necessary.

5. Disposal: Use any included waste bags to dispose of waste material safely and efficiently

